

3.2 Technical data

Drive pressure	2.5 (29) to 4.2 (44) bar (psi)
Air consumption	20 – 40 NL/min
Water use	30 - 50 ml/min 3.2.3
Frequency	6 – 6.5 kHz 3.2.2
Recommended pressure	0.1 - 2 N

The SONICflex can be mounted on all MULTIflex couplings. [3.2.7](#)

- Oscillation level 1 = 120 +/- 15 µm
- Oscillation level 2 = 160 +/- 15 µm
- Oscillation level 3 = 240 +/- 15 µm